



## EEG BOARD ADAPTER CODES

<b>WHITE TIPS</b>	<b>WIRE COLOR</b>	<b>RED TIPS</b>
FP1	<b>BROWN</b>	FP2
F3	<b>RED</b>	F4
C3	<b>ORANGE</b>	C4
P3	<b>YELLOW</b>	P4
O1	<b>GREEN</b>	O2
F7	<b>BLUE</b>	F8
T3	<b>PURPLE</b>	T4
T5	<b>GRAY</b>	T6
GND	<b>WHITE</b>	CZ
FZ	<b>BLACK</b>	PZ